


# For Rapid Airway Placement

The McMurray Enhanced Airway (MEA)

OPEN • PLACE • GO



In emergency response, speed and accuracy count to reduce hypoxia.



## Intuitive, fast and easy to place

The MEA quickly opens obstructed airways. The MEA's small, flexible tubing can accommodate narrow airway openings. Simply place the MEA in the midline of the mouth, follow the hard palate, and slide between the molars.

## Excellent intraoral ventilation option

Mask ventilation can be challenging, especially with obese, edentulous or bearded patients. Intraoral ventilation helps to bridge prior to securing an airway. Place the MEA and use the optional connector to connect to a resuscitator bag, then ventilate while manually closing the nose and mouth. The MEA is an excellent option for quickly establishing a patent airway.

## Provides supplemental oxygen during intubation

To decrease hypoxia during intubation, place MEA on the left side of the mouth, attach oxygen source, and intubate around it.

For more information or to place an order visit:  
[www.mcmurraymed.com](http://www.mcmurraymed.com)

Roxanne McMurray,  
Anesthesia Provider and  
Inventor of the MEA

“The MEA is designed to be fast and easy to use and to maintain a patent airway for a wide variety of patients.”

# The MEA simplifies airway management for emergency medical responders.

The MEA is a distal pharyngeal airway (DPA) designed to open and maintain a patient's upper airway beyond the base of the tongue.

## Soft, flexible material and longer tubing than other currently available oral airways




The MEA's flexible material is intended to ease insertion, helping reduce coughing, gagging and oral/dental damage. Longer tubing helps displace the redundant pharyngeal tissue that often blocks sedated or unconscious patients' airways. This helps establish and maintain a patent airway for improved oxygenation and ventilation.

## Eliminates nose bleeds due to nasal airway insertion

Anti-coagulated patients or those on an aspirin regimen are at risk of nosebleeds if a nasal airway is placed. The MEA is placed orally, avoiding this risk.

## Compact, individually packaged

Individually packaged in an easy-open poly bag. Single-patient use. Sized to accommodate a broad adult patient population.



## Contact Us

Visit [mcmurraymed.com](http://mcmurraymed.com) for more information, to see videos of the MEA in use, or to request a sample.

email: [info@mcmurraymed.com](mailto:info@mcmurraymed.com)

